

## Maintenance of *Daphnia pulex* Cultures

### STEPS

### COMMENTS

#### Scope and Application

This method is based on EPA 821-R-02-012 *Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms*, 5<sup>th</sup> Ed. October 2002 and the test organism culturing experience of laboratory personnel. It is used to ensure that *Daphnia pulex* bioassay test organism cultures are maintained so as to provide suitable test organisms should the laboratory begin testing with *D. pulex*.

#### Summary of the Method

Cultured organisms are fed and transferred on a regular basis to prepared vessels containing fresh media. Health of the organisms is monitored and recorded in the *Culture Records for C. dubia, D. pulex and H. azteca* logbook.

#### Interferences

Improperly prepared culture vessels, culture media and food can be detrimental to the survival and reproduction of cultured *D. pulex*.

#### Materials

##### 1. *Daphnia pulex* “Mass” culture

*Daphnia pulex* are not maintained as a source of test organisms so a “test species” culture record is not maintained at this time; however, culture feeding and maintenance is recorded in the “Culture records for...” laboratory notebook. The “Mass” culture is kept in a 4 L beaker on the culture tray labeled “B4”. The *Daphnia pulex* “Reserve” culture is kept in an 8” (20.32 cm) Carolina<sup>®</sup> dish on Tray B4.

##### 2. *Daphnia pulex* “Reserve” culture

##### 3. Well water

##### 4. Siphon tube

##### 5. One 4-L beaker

##### 6. Carolina dish, 8” (20.32 cm)

##### 7. Wastewater bucket, 5-gal (18.93 L)

##### 8. Clean large bore transfer pipette

##### 9. *Culture Records for C. dubia, D. pulex and H. azteca* laboratory notebook.

#### Methods

##### *Mass Cultures*

##### 1. Remove the “mass culture” from tray “B4” and place it on the bioassay work table. Siphon off approximately half of the water and *D. pulex* from the “current” culture into the waste bucket.

Clean once per week (Friday) unless culture conditions warrant cleaning more often (algae blooms, culture crash, cloudiness, etc.).

##### 2. Pour the remaining culture water and *D. pulex* into the clean 4 L beaker.

##### 3. Label the new, clean beaker with “*Daphnia pulex* Mass Culture”, the “date”, and your “initials”.

STEPS	COMMENTS
<ol style="list-style-type: none"> <li>4. Add well water to the new beaker to bring the volume to the 3-L mark. Return the new culture to tray “B4”.</li> <li>5. Feed the culture per SOP TA-02.00.</li> <li>6. Dispose of the organisms in the waste bucket per Step 7 of the “Reserve Culture” method.</li> <li>7. Place the old culture beaker on the dirty glassware cart for cleaning.</li> <li>8. Record culture maintenance and feeding in the <i>Culture Records for C. dubia, D. pulex and H. azteca</i> laboratory notebook.</li> </ol>	See page 3
<b><i>Reserve Cultures</i></b>	
<ol style="list-style-type: none"> <li>1. Fill the clean 8” (20.32 cm) Carolina dish 3/4 full with well water.</li> <li>2. Remove the “reserve culture” from tray “B4” and place it on the bioassay work table. Using the transfer pipette, carefully transfer 50-75 of the <i>D. pulex</i> to the new, clean dish.</li> <li>3. Label the new, clean dish with “<i>Daphnia pulex</i> Reserve Culture”, the “date”, and your “initials”.</li> <li>4. Place the new “reserve culture” on tray “B4”.</li> <li>5. Feed the culture per SOP TA-02.00.</li> <li>6. Pour the culture water and remaining daphnia in the “old” “reserve dish” into the waste bucket.</li> <li>7. Euthanize the organisms in the waste bucket by filling the bucket with ice and letting it sit for 20-30 minutes. Pour the iced water and euthanized organisms into the laboratory sink and flush the sink with tap water.</li> <li>8. Place the old culture dish on the dirty glassware cart for cleaning.</li> <li>9. Record culture maintenance and feeding in the <i>Culture Records for C. dubia, D. pulex and H. azteca</i> laboratory notebook.</li> </ol>	<p>Clean once per week (Friday) unless culture conditions warrant cleaning more often (algae blooms, culture crash, cloudiness, etc.).</p> <p>Transfer both large and small organisms (approximately half large and half small).</p>

## **STEPS**

## **COMMENTS**

### **References**

1. EPA 821-R-02-012 *Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms*, 5<sup>th</sup> Ed. October 2002.
2. Feeding Test-Organism Cultures, DEP SOP TA-02.00
3. Biology Quality Manual, FDEP Bureau of Laboratories, March 2008

### **Appendix of Changes**

9/25/2008 Added *Scope and Application, Summary of the Method, Interferences, References and Appendix of Changes*